



CALGIDRESS™

Sterile Calcium Alginate Wound Dressing

High Absorbency - Accelerates Wound Healing

INTRODUCTION

Calcium alginate dressings are primarily used for managing moderately to heavily exuding wounds, helping to maintain a moist wound environment and promote healing. They are particularly useful for wounds with irregular shapes and significant tissue damage. They can also be used for Diabetic Foot ulcers pressure ulcers, vascular ulcers, donor sites, and surgical incisions.

MODE OF ACTION

- ✓ **PROMOTS MOIST WOUND HEALING**
- ✓ **MANAGING EXUDATE**
- ✓ **CONTROLS BLEEDING**
- ✓ **FACILITATES TISSUE REGENERATION**
- ✓ **SUPPORT FOR WOUND CLOSURE**
- ✓ **EASE OF REMOVAL**

