

HAEMOSPON™

Sterile Haemostatic Absorbable Gelatin Sponge

PRECAUTIONS & WARNINGS

- **HAEMOSPON** is a sterile product and should not be re-sterilized.
- Sterility not guaranteed if the package is damaged or tampered.
- Discard the unused portions of sponge.
- Use of **HAEMOSPON** is not recommended for infected sites.
- By absorbing the fluid, the sponge may expand and thus impinge on neighboring structures. Hence, careful decompression is advised when packaging into cavities.
- While using **HAEMOSPON**, other normal haemostatic procedures including suture, pressure and ligation should be employed.
- Avoid in patients hypersensitive to bovine materials.

STORAGE

- **HAEMOSPON** must be stored in a dry place with a shelf life of 3 years from the date of manufacture and the expiry date is printed on the packaging.



HAEMOSPON Size: 80 x 40 x 10 mm Pack of 5 Units



KRATI
BIOTECH
Healing People

info@kratibiotech.com
kratibiotech2009@gmail.com
@ 894302880 / 805721751
website : www.kratibiotech.com